

immu

THE SOLUTION

Dietary Supplement

"Revolutionary micro-sized emulsion that supports a healthy immune system"



How to Use:

Shake well to homogenize the emulsion (product separation is typical and no quality defect).

- Adult Serving Size: 0.2 ml (4 drops) daily. (maximum dose 0.5 ml daily)
- Child Serving Size: 0.1 ml (2 drops) once a week. (maximum dose 0.25 ml daily)

Ingredients (per 2 ml):

Chondroitin Sulfate & Phosphatidylcholine (13mg), Vitamin D as Cholecalciferol (10µg), Water, Salt.

"The secret is in the molecular bonding"

0.1 fl oz (3 ml)

imno

THE SOLUTION

How to Load Needleless Syringe

- Prise the plastic flip-top disk from the top of the vial.
- Remove the vial adapter from the packaging and push it onto the top of the vial. The plastic spike will pierce the grey stopper visible in the centre of the metal top, and the adapter will clip into place.
- Push the closed syringe firmly into the top of the vial adapter (depressing the blue rubber valve).
- Invert the vial with syringe underneath and draw the syringe open to the required amount (i.e. 0.2ml).
- Turn the vial upright and remove the syringe from the vial adapter. The valve will close and seal the vial/adapter on syringe removal.
- Use the required dose from the needleless syringe.

Note: Every division on the supplied 1ml needleless syringe equals 0.02ml, with major divisions numbered at each 0.1ml.



Transport: room temperature | Protect: from light
Store: refrigerated.

Designed by Molecular Biologist, Professor Marco Ruggiero, MD. PhD.

Made in New Zealand for imuno Corporation
Lini Highway, Port Vila, Vanuatu.